

## EBLV Ab ELISA 480

**EAN Code:** 8595635302596

**Catalog number:** EBL480

**Package size:** 480 tests

**Storage:** 2-8 °C

**Producer:** TestLine Clinical Diagnostics s.r.o.



### Description:

- 2 versions of the kit (screening, confirmative).
- Microtitre wells are coated with the virus antigen (p-24) using monoclonal antibodies (screening microtitre plate is coated with the EBLV antigen).
- The proportion of screening and confirmative microplates can be changed on request.
- Screening kit (5 microplates – 12 × 8 wells) for up to 480 screening tests for samples of blood serum, including controls.
- The whole kit may be used sequentially for smaller batches of samples.
- The positive control corresponds to the Danish positive reference serum E4 (diluted 10 times).
- Total time of survey according to sample incubation duration:
  - over night: results in second day.
  - rapid test: approximately 2.5 h.
- The evaluation: comparing the absorbance of the sample against the absorbance of the negative control.
- Tested by the Veterinary Research Institute, Brno, Czech Republic and the State Veterinary Institutes in the Czech Republic.
- Approved by the Institute for State Control of Veterinary Biologicals and Medicaments, Czech Republic.
- The kit is manufactured in accordance with Good Manufacturing Practice (GMP).
- Shelf life: 12 months.

### Applications:

- Screening of EBL occurrence in a cattle population.

- Control of current sanitation process in cattle herds.
- Ongoing screening in herds after sanitation.
- Control of transported and quarantined animals (export, import).

**Brief assay procedure:**

1. Dilute samples and controls (1:20) in microplate wells.
2. Incubate:
  - a) over night: 14–18 hours at 4–8 °C
  - b) rapid test: 60 min at 37 °C
3. Aspirate and wash the wells 4 times.
4. Add Conjugate.
5. Incubate for 60 min at 37 °C.
6. Aspirate and wash the wells 4 times.
7. Add substrate (TMB-Complete).
8. Incubate for 10 min at room temperature.
9. Add Stopping solution (H<sub>2</sub>SO<sub>4</sub>).
10. Read the absorbance photometrically at 450 nm.
11. Evaluate results.