

## BHV-1 Ab ELISA 480

**EAN Code:** 8595635302572

**Catalog number:** BHd480

**Package size:** 480 tests

**Storage:** 2-8 °C

**Producer:** TestLine Clinical Diagnostics s.r.o.



### Description:

- The kit is standardized according to the OIE International Standard Serum EU2.
- The kit enables 480 screening tests, including controls.
- The whole kit may be used at one time or sequentially for smaller batches of samples.
- Microtitre wells are coated with the purified and inactivated BHV-1 antigen.
- Total assay time: approximately 2 h.
- The evaluation is based on ratio of the absorbance of the sample to the absorbance of the positive control serum-CUT OFF expressed in percentage (%).
- Tested by the Veterinary Research Institute, Brno, Czech Republic and State Veterinary Institutes in the Czech Republic.
- Approved by the Institute for State Control of Veterinary Biologicals and Medicaments, Czech Republic.
- The kit is manufactured in accordance with Good Manufacturing Practice (GMP).
- Shelf life: 15 months.

### Application:

- Screening of IBR occurrence in a cattle population.
- Identification of infected or vaccinated animals, checking of vaccination effectiveness.
- Obtaining the background information to formulate sanitation programmes for cattle herds or regions, control of current sanitation process in cattle herds.
- Evaluation of effects of the sanitation programme, ongoing screening in herds after sanitation.
- Control of transported and quarantined animals (export, import).

### Brief assay procedure:

1. Dilute samples and controls (1:20) in microplate wells.
2. Incubate for 60 min at 37 °C.
3. Aspirate and wash the wells 4 times.
4. Add Conjugate.
5. Incubate for 30 min at 37 °C.
6. Aspirate and wash the wells 4 times.
7. Add substrate (TMB-Complete).
8. Incubate for 15 min at room temperature.
9. Add Stopping solution (H<sub>2</sub>SO<sub>4</sub>).
10. Read absorbance photometrically at 450 nm.
11. Evaluate results.