

BHV-1 Ab ELISA 480



EAN Code: 8595635302572 Catalog number: BHd480 Package size: 480 tests

Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.

Description:

- The kit is standardized according to the OIE International Standard Serum EU2.
- The kit enables 480 screening tests, including controls.
- The whole kit may be used at one time or sequentially for smaller batches of samples.
- Microtitre wells are coated with the purified and inactivated BHV-1 antigen.
- Total assay time: approximately 2 h.
- The evaluation is based on ratio of the absorbance of the sample to the absorbance of the positive control serum-CUT OFF expressed in percentage (%).
- Tested by the Veterinary Research Institute, Brno, Czech Republic and State Veterinary Institutes in the Czech Republic.
- Approved by the Institute for State Control of Veterinary Biologicals and Medicaments, Czech Republic.
- The kit is manufactured in accordance with Good Manufacturing Practice (GMP).
- Shelf life: 15 months.

Application:

- Screening of IBR occurrence in a cattle population.
- Identification of infected or vaccinated animals, checking of vaccination effectiveness.
- Obtaining the background information to formulate sanitation programmes for cattle herds or regions, control of current sanitation process in cattle herds.
- Evaluation of effects of the sanitation programme, ongoing screening in herds after sanitation.
- Control of transported and quarantined animals (export, import).

Brief assay procedure:

Exportováno 22. 11. 2025 Stránka 1 / 2



- 1. Dilute samples and controls (1:20) in microplate wells.
- 2. Incubate for 60 min at 37 °C.
- 3. Aspirate and wash the wells 4 times.
- 4. Add Conjugate.
- 5. Incubate for 30 min at 37 °C.
- 6. Aspirate and wash the wells 4 times.
- 7. Add substrate (TMB-Complete).
- 8. Incubate for 15 min at room temperature.
- 9. Add Stopping solution (H2SO4).
- 10. Read absorbance photometrically at 450 nm.
- 11. Evaluate results.

Exportováno 22. 11. 2025 Stránka 2 / 2