

BlueBLOT-LINE Helicobacter IgG



EAN Code: 8595635307805

Catalog number: BD-HpGL24

Package size: 24 tests

Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.

Description

- Each single test contains all the reagents (ready-to-use) necessary for performing the test
- Antigen used: combination of selected parts of the specific antigens of *H. pylori* (CagA, VacA, UreA, NAP, HpaA, HcpC, GroEL)
- Recombinant antigens are transferred to the nitrocellulose membrane, if specific antibodies are present in the sample, they will bind to the respective antigens – the complex is labelled with Conjugate and detected through a color reaction with substrate
- Short incubation periods, total assay time: approximately 1.25 h.
- Evaluation using Immunoblot Software (TestLine)

Applications

- Diagnostics of *H. pylori*
- Confirmation for ELISA tests

Brief assay procedure:

1. Insert testing strips and cartridges into the BlueDiver Instrument.
2. Pipette samples.
3. Fully automated incubation and washing steps in BlueDiver Instrument
4. Incubation 30 minutes.
5. Wash 3 x 2 min.
6. Conjugate dispensing.
7. Incubation 20 min.
8. Wash 3 x 2 min.
9. Substrate dispensing.
10. Incubation 10 min.
11. Wash 1x 2 min.
12. Evaluation of test results using the BlueScan scanner and Immunoblot Software.