

EBLV Ab ELISA 480

EAN Code: 8595635302596

Catalog number: EBL480

Package size: 480 tests

Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.



Description:

- Microtitre wells are coated with the virus antigen (p-24) using monoclonal antibodies (screening microtitre plate is coated with the EBLV antigen).
- The proportion of screening and confirmative microplates can be changed on request.
- Screening kit (5 microplates – 12 × 8 wells) for up to 480 screening tests for samples of blood serum, including controls.
- The whole kit may be used sequentially for smaller batches of samples.
- The positive control corresponds to the Danish positive reference serum E05 (diluted 10 times).
- Total time of survey according to sample incubation duration:
 - over night: results in second day.
 - rapid test: approximately 2.5 h.
- The evaluation: comparing the absorbance of the sample against the absorbance of the negative control.
- Tested by the Veterinary Research Institute, Brno, Czech Republic and the State Veterinary Institutes in the Czech Republic.
- Approved by the Institute for State Control of Veterinary Biologicals and Medicaments, Czech Republic.
- The kit is manufactured in accordance with Good Manufacturing Practice (GMP).
- Shelf life: 12 months.

Applications:

- Screening of EBL occurrence in a cattle population.
- Control of current sanitation process in cattle herds.

- Ongoing screening in herds after sanitation.
- Control of transported and quarantined animals (export, import).

Brief assay procedure:

1. Dilute samples and controls (1:20) in microplate wells.
2. Incubate:
 - a) over night: 14–18 hours at 4–8 °C
 - b) rapid test: 60 min at 37 °C
3. Aspirate and wash the wells 4 times.
4. Add Conjugate.
5. Incubate for 60 min at 37 °C.
6. Aspirate and wash the wells 4 times.
7. Add substrate (TMB-Complete).
8. Incubate for 10 min at room temperature.
9. Add Stopping solution (H₂SO₄).
10. Read the absorbance photometrically at 450 nm.
11. Evaluate results.