


Determination of antibodies against SARS-CoV-2

Correlation of VNT and ELISA kit results
Compliance with WHO standard

A standard range of methodological approaches is available for the detection of antibodies against SARS-CoV-2 antigens. For the detection of antibodies, commercial tests are available to detect IgA, IgG, and IgM based on immuno-enzymatic principles such as ELISA, Immunoblot, CLIA.

The antigen in the assays is usually either a structural nucleocapsid (NP) protein or a spike (S) protein (or its S1 part or only the Receptor-Binding Domain – RBD).



According to the latest literature, it is most suitable to use a virus neutralization test (VNT) to determine the protective activity of antibodies present in a patient's serum. However, the determination of neutralizing antibodies (NAbs) by the VNT method is time-consuming (2-4 days) and, in addition, requires work with a live virus (BSL3 protection level laboratory). Standard commercially available assays detect anti-SARS-Cov-2 binding antibodies (BAbs) and do not have the ability to differentiate between NAbs and BAbs.

The correlation of the results between VNT and ELISA was performed in a comparative study of 100 samples.

The TestLine EIA COVID-19 RBD diagnostic kits were used for testing, where a significant degree of agreement was found. Both semi-quantitative evaluation of IgG results in the Positivity Index (IP) and the quantitative evaluation corresponds to the VNT titer.

| VNT titer | Elisa – TESTLINE RBD | | | | | | ELISA – competitor DE | | | VNT titer | Elisa – TESTLINE RBD | | | | | | ELISA – competitor DE | | |
|-----------|----------------------|--------|------|-------|------|--------|-----------------------|------|------|-----------|----------------------|--------|------|-------|------|--------|-----------------------|------|------|
| | IgG | | IgM | | IgA | | IgG | IgM | IgA | | IgG | | IgM | | IgA | | IgG | IgM | IgA |
| | IP | U/ml | IP | U/ml | IP | U/ml | IP | IP | IP | IP | U/ml | IP | U/ml | IP | U/ml | IP | IP | IP | |
| 20 | 3.83 | 97.3 | 0.67 | 15 | 0.72 | 5.57 | 1.29 | neg. | 0.93 | 80 | 2.50 | 59.81 | 0.63 | 14.39 | 3.92 | 204.5 | 1.35 | 1.76 | 1.95 |
| 20 | 0.48 | 11.02 | 0.35 | 10.04 | 0.55 | 3.03 | neg. | neg. | neg. | 80 | 3.40 | 84.14 | 0.72 | 15.78 | 1.60 | 50.83 | 1.16 | 0.90 | 2.61 |
| 20 | 2.82 | 68.18 | 0.45 | 11.56 | 0.78 | 6.55 | neg. | neg. | 0.94 | 80 | 8.47 | 238.15 | 0.63 | 14.34 | 3.54 | 177.03 | 1.96 | neg. | 3.10 |
| 20 | 2.57 | 61.7 | 0.54 | 12.91 | 0.35 | 9.82 | neg. | neg. | neg. | 80 | 3.28 | 80.38 | 0.56 | 13.2 | 1.60 | 50.45 | 1.28 | neg. | 1.47 |
| 20 | 1.26 | 26.81 | 0.64 | 14.51 | 1.67 | 54.16 | neg. | neg. | 0.94 | 80 | 8.07 | 226.12 | 1.92 | 63.72 | 9.44 | 320 | 1.46 | 1.48 | 8.97 |
| 20 | 1.48 | 47.1 | 0.34 | 9.99 | 1.05 | 21.56 | 1.04 | neg. | 0.97 | 80 | 2.88 | 69.9 | 0.85 | 17.75 | 9.44 | 320 | 1.01 | neg. | 2.45 |
| 20 | 2.45 | 117.27 | 1.59 | 47.85 | 1.33 | 38.22 | 1.35 | neg. | 2.84 | 80 | 3.40 | 84.23 | 0.64 | 14.43 | 3.41 | 167.58 | 1.28 | 2.08 | 2.90 |
| 40 | 2.58 | 61.86 | 0.33 | 9.8 | 0.85 | 7.73 | 1.20 | neg. | 1.09 | 80 | 9.78 | 277.83 | 0.87 | 17.99 | 1.05 | 21.3 | 2.31 | neg. | 5.53 |
| 40 | 1.48 | 32.8 | 0.22 | 8.11 | 0.69 | 5.18 | 1.35 | neg. | neg. | 80 | 8.49 | 238.71 | 1.09 | 24.05 | 2.26 | 85.71 | 1.54 | 0.94 | 1.43 |
| 40 | 4.96 | 131.53 | 1.30 | 34.45 | 0.92 | 18.83 | 2.71 | neg. | 1.05 | 80 | 3.72 | 93.82 | 1.03 | 21.27 | 0.46 | 1.62 | 1.01 | 2.72 | neg. |
| 40 | 4.11 | 105.57 | 0.43 | 11.23 | 0.39 | 0.4 | 1.51 | neg. | neg. | 80 | 10.31 | 293.91 | 0.75 | 16.19 | 6.66 | 320 | 2.37 | neg. | 1.91 |
| 40 | 4.11 | 105.67 | 0.22 | 7.99 | 1.49 | 44.95 | 1.28 | neg. | 1.10 | 80 | 9.44 | 267.67 | 0.44 | 11.43 | 1.54 | 47.76 | 1.75 | neg. | 1.61 |
| 40 | 4.79 | 126.35 | 0.49 | 12.21 | 1.94 | 67.97 | 0.99 | neg. | 1.00 | 80 | 3.40 | 208.56 | 1.24 | 31.33 | 1.62 | 54 | 1.94 | neg. | 1.38 |
| 40 | 7.86 | 219.53 | 0.56 | 13.28 | 5.80 | 320 | 2.72 | neg. | 3.60 | 80 | 4.22 | 287.25 | 0.43 | 11.36 | 0.64 | 14.48 | 2.24 | neg. | 1.00 |
| 40 | 1.84 | 42.33 | 0.56 | 13.24 | 0.49 | 2.01 | neg. | neg. | neg. | 80 | 4.62 | 320 | 0.66 | 14.82 | 1.32 | 37.79 | 3.76 | 0.99 | 1.65 |
| 40 | 8.74 | 246.24 | 0.39 | 10.61 | 3.11 | 146.71 | 2.35 | neg. | 1.92 | 80 | 4.33 | 297.59 | 0.27 | 8.92 | 0.78 | 16.65 | 3.09 | 2.03 | 1.09 |
| 40 | 3.39 | 83.76 | 0.52 | 12.66 | 1.34 | 37.4 | 1.47 | neg. | neg. | 160 | 11.94 | 320 | 1.61 | 49.02 | 3.14 | 148.85 | 6.16 | 1.98 | 2.72 |
| 40 | 4.03 | 103.22 | 0.26 | 8.73 | 1.59 | 50.19 | 1.43 | 1.98 | - | 160 | 11.94 | 320 | 0.88 | 18.11 | 1.88 | 65.03 | 6.38 | - | - |
| 40 | 3.34 | 82.44 | 0.52 | 12.62 | 0.68 | 5.03 | neg. | neg. | neg. | 160 | 11.42 | 320 | 0.21 | 7.83 | 2.44 | 99.09 | 5.12 | - | - |
| 40 | 2.57 | 61.61 | 0.50 | 12.34 | 1.84 | 62.73 | neg. | 1.27 | 0.92 | 160 | 9.43 | 267.3 | 1.48 | 42.56 | 1.87 | 64.39 | 3.81 | 1.79 | 2.07 |
| 40 | 3.81 | 96.74 | 0.51 | 12.5 | 0.84 | 7.53 | 1.47 | neg. | 0.99 | 160 | 5.42 | 145.54 | 1.19 | 29.12 | 2.53 | 104.97 | 1.74 | neg. | 1.37 |
| 40 | 6.17 | 168.38 | 0.32 | 9.51 | 1.45 | 43.16 | 1.80 | neg. | 1.78 | 160 | 9.83 | 279.43 | 0.32 | 9.59 | 2.57 | 108.18 | 3.02 | 1.08 | 3.26 |
| 40 | 2.39 | 56.77 | 1.64 | 50.41 | 0.83 | 7.42 | 0.95 | neg. | neg. | 160 | 12.06 | 320 | 1.22 | 30.64 | 6.51 | 320 | 5.02 | 0.92 | 7.55 |
| 40 | 4.57 | 320 | 0.85 | 17.72 | 2.59 | 125.8 | 2.63 | 1.41 | 2.10 | 160 | 9.86 | 280.46 | 0.61 | 14.1 | 1.51 | 45.97 | 3.21 | neg. | 3.59 |
| 40 | 3.82 | 248.14 | 1.11 | 25.07 | 0.72 | 15.67 | 1.45 | neg. | 1.60 | 160 | 11.76 | 320 | 0.54 | 12.95 | 3.29 | 159.2 | 3.56 | neg. | 2.69 |
| 80 | 11.25 | 320 | 0.64 | 14.55 | 1.50 | 45.59 | 4.90 | - | - | 160 | 11.55 | 320 | 0.53 | 12.75 | 1.37 | 38.68 | 3.64 | neg. | 1.54 |
| 80 | 11.38 | 320 | 0.86 | 17.87 | 4.18 | 222.34 | 4.13 | - | - | 160 | 10.23 | 291.65 | 0.68 | 15.16 | 1.09 | 21.61 | 3.35 | neg. | 1.69 |
| 80 | 6.19 | 168.76 | 2.25 | 79.43 | 0.67 | 4.91 | 2.12 | neg. | 0.98 | 160 | 11.41 | 320 | 0.44 | 11.48 | 1.36 | 38.29 | 3.90 | neg. | 2.50 |
| 80 | 11.79 | 320 | 0.88 | 18.16 | 2.35 | 92.49 | 5.43 | neg. | 3.99 | 160 | 11.08 | 317.32 | 0.92 | 18.81 | 2.18 | 80.36 | 2.39 | neg. | 2.07 |
| 80 | 4.23 | 109.43 | 0.15 | 7.01 | 3.20 | 152.78 | 1.92 | neg. | 7.23 | 160 | 2.23 | 52.59 | 1.22 | 30.52 | 5.98 | 320 | 0.95 | neg. | 3.62 |
| 80 | 5.78 | 156.44 | 0.47 | 11.93 | 0.76 | 6.2 | 2.56 | neg. | 3.38 | 160 | 9.11 | 257.61 | 0.78 | 16.68 | 1.38 | 39.45 | 2.16 | neg. | 1.01 |
| 80 | 6.69 | 184.18 | 0.92 | 18.73 | 2.78 | 122.99 | 2.73 | neg. | 2.47 | 160 | 12.01 | 320 | 1.06 | 22.79 | 4.17 | 221.63 | 2.68 | 2.89 | 4.97 |
| 80 | 4.91 | 129.93 | 0.45 | 11.6 | 1.04 | 18.79 | 1.39 | neg. | 1.03 | 160 | 12.04 | 320 | 0.81 | 17.09 | 1.80 | 60.68 | 4.42 | neg. | 1.93 |
| 80 | 2.80 | 67.61 | 0.36 | 10.2 | 1.93 | 67.72 | 1.35 | neg. | neg. | 160 | 11.55 | 320 | 2.40 | 90.06 | 1.44 | 42.39 | 2.39 | neg. | 0.95 |
| 80 | 5.02 | 133.5 | 0.35 | 10.04 | 0.33 | 9.54 | 2.41 | neg. | neg. | 160 | 11.69 | 320 | 0.93 | 18.93 | 3.36 | 320 | 4.34 | 3.43 | 2.47 |
| 80 | 4.63 | 121.56 | 0.33 | 9.75 | 1.14 | 27.16 | 1.71 | neg. | 0.99 | 160 | 10.83 | 309.89 | 0.68 | 15.08 | 2.57 | 107.83 | 1.61 | neg. | 1.90 |

| VNT titer | Elisa – TESTLINE RBD | | | | | | ELISA – competitor DE | | | VNT titer | Elisa – TESTLINE RBD | | | | | | ELISA – competitor DE | | | |
|-----------|----------------------|--------|------|--------|------|--------|-----------------------|------|------|-----------|----------------------|--------|--------|-------|--------|--------|-----------------------|------|------|------|
| | IgG | | IgM | | IgA | | IgG | IgM | IgA | | IgG | | IgM | | IgA | | IgG | IgM | IgA | |
| | IP | U/ml | IP | U/ml | IP | U/ml | IP | IP | IP | IP | U/ml | IP | U/ml | IP | U/ml | IP | IP | IP | | |
| 160 | 11.81 | 320 | 2.18 | 76.01 | 4.52 | 246.6 | 5.12 | 1.17 | 3.06 | 320 | 12.23 | 320 | 1.43 | 40.4 | 9.26 | 320 | 3.98 | 1.12 | 4.81 | |
| 160 | 12.01 | 320 | 4.16 | 227.17 | 9.23 | 320 | 1.97 | 0.99 | 8.45 | 320 | 11.01 | 315.25 | 6.40 | 4.16 | 221.45 | 3.69 | 3.80 | 2.57 | | |
| 160 | 3.70 | 93.35 | 5.75 | 320 | 7.80 | 320 | neg. | 3.51 | 3.69 | 320 | 11.42 | 320 | 0.85 | 17.7 | 7.17 | 320 | 3.89 | neg. | 7.47 | |
| 160 | 4.81 | 320 | 7.22 | 320 | 3.65 | 222.5 | 3.63 | neg. | 5.07 | 320 | 11.69 | 320 | 0.45 | 11.56 | 2.66 | 114.25 | 2.41 | neg. | 2.11 | |
| 160 | 4.68 | 320 | 0.37 | 10.41 | 1.72 | 59.44 | 4.14 | neg. | 3.94 | 320 | 12.06 | 320 | 1.83 | 59.16 | 5.30 | 302.25 | 5.72 | 2.01 | 5.24 | |
| 160 | 4.46 | 309.84 | 2.03 | 68.61 | 2.51 | 118.1 | 3.86 | 0.98 | 5.41 | 320 | 11.69 | 320 | 7.06 | 320 | 9.32 | 320 | 5.53 | 4.70 | 4.40 | |
| 160 | 4.62 | 320 | 0.42 | 11.09 | 2.30 | 98.91 | 5.13 | neg. | 4.42 | 320 | 12.23 | 320 | 0.78 | 16.56 | 3.63 | 183.28 | 6.52 | neg. | 4.84 | |
| 320 | 12.01 | 320 | 3.62 | 185 | 8.88 | 320 | 7.52 | - | - | 320 | 11.81 | 320 | 1.91 | 63.08 | 2.69 | 116.39 | 2.83 | 3.32 | 2.04 | |
| 320 | 5.80 | 157.62 | 1.30 | 34.45 | 8.93 | 320 | 2.74 | 1.21 | 5.54 | 320 | 640 | 11.94 | 320 | 0.99 | 19.88 | 8.73 | 320 | 6.93 | neg. | 5.31 |
| 320 | 5.20 | 138.86 | 1.62 | 49.52 | 1.38 | 39.57 | 1.49 | neg. | 1.35 | 320 | 640 | 11.76 | 320 | 1.25 | 31.66 | 7.65 | 320 | 6.16 | neg. | 7.66 |
| 320 | 11.49 | 320 | 0.98 | 19.63 | 3.65 | 184.7 | 5.29 | neg. | 5.17 | 320 | 640 | 10.08 | 286.95 | 2.72 | 114.9 | 0.32 | 9.35 | 4.26 | 2.14 | 0.98 |
| 320 | 12.23 | 320 | 0.90 | 18.4 | 8.72 | 320 | 4.28 | neg. | 2.34 | 320 | 640 | 11.94 | 320 | 2.25 | 79.18 | 5.79 | 320 | 5.40 | 2.10 | 4.20 |
| 320 | 11.88 | 320 | 0.29 | 9.14 | 3.04 | 141.72 | 4.12 | neg. | 3.78 | 320 | 640 | 12.01 | 320 | 9.57 | 320 | 4.74 | 262.3 | 4.29 | 3.09 | 4.02 |
| 320 | 11.88 | 320 | 0.67 | 15 | 2.24 | 84.64 | 5.08 | 1.10 | 3.54 | 320 | 640 | 12.01 | 320 | 9.57 | 320 | 4.74 | 262.3 | 4.29 | 3.09 | 4.02 |
| 320 | 11.88 | 320 | 0.67 | 15 | 2.24 | 84.64 | 5.08 | 1.10 | 3.54 | 320 | 1280 | 11.40 | 320 | 7.93 | 320 | 9.39 | 320 | 5.36 | 3.30 | 7.51 |

Results interpretation

TL IP: light blue negative: < 0,9 medium blue borderline: 0,9 - 1,1 dark blue positive: > 1,1
 TL U/ml: light blue negative: < 18 medium blue borderline: 18 - 22 dark blue positive: > 22
 competitor DE: light blue negative: < 0,8 medium blue borderline: 0,8 - 1,1 dark blue positive: > 1,1

The correlation of results with VNT method

VNT vs EIA TESTLINE IgG

| VNT | EIA | |
|------------------|-------------|-----|
| | pos | neg |
| | pos | 99 |
| neg | 0 | 0 |
| Agreement | 99 % | |

All classes of VNT antibodies vs EIA TESTLINE

| VNT | EIA | |
|------------------|-------------|-----|
| | pos | neg |
| | pos | 99 |
| neg | 0 | 0 |
| Agreement | 99 % | |


VNT vs EIA competitor IgG (DE)

| VNT | EIA | |
|------------------|-------------|-----|
| | pos | neg |
| | pos | 92 |
| neg | 0 | 0 |
| Agreement | 92 % | |


All classes of VNT antibodies vs EIA competitor (DE)

| VNT | EIA | |
|------------------|-------------|-----|
| | pos | neg |
| | pos | 96 |
| neg | 0 | 0 |
| Agreement | 96 % | |

Mean Positive Index (IP) values of IgG anti-RBD antibodies (TestLine) in relation to individual VNT titers




Mean values of U/ml IgG anti-RBD antibodies (TestLine) in relation to individual VNT titers



Compliance with WHO standard


Titration of the WHO 20/136 IgG standard

A group of convalescent plasma from patients recovered from COVID-19, containing high titers of antibodies against SARS-CoV-2.




Titration of the WHO 20/136 IgA standard

A group of convalescent plasma from patients recovered from COVID-19, containing high titers of antibodies against SARS-CoV-2.




Titration of the NIBSC 20/162 IgG standard

High titer anti-SARS-CoV-2 antibody material was used to assess and compare relative sensitivities for Anti-Sars CoV-2 determination of antibodies by dilution at the endpoint.



Titration of the NIBSC 20/162 IgA standard

High titer anti-SARS-CoV-2 antibody material was used to assess and compare relative sensitivities for Anti-Sars CoV-2 determination of antibodies by dilution at the endpoint.



The WHO 20/136 International Standard for Anti-SARS-CoV-2 is intended for the calibration and harmonization of serological tests detecting the anti-SARS-CoV-2 neutralizing effect and for determining the validated level of binding antibodies. The NIBSC Anti-SARS-Cov-2 standard is designated as Quality Control Reagent and monitors the sensitivity of serological tests.

By titrating WHO 20/136 Standard, TestLine responded to the requirement to accurately determine the level of antibodies after vaccination. According to the achieved titer of the Standard was recommended further dilution to determine the level of antibody decline in strongly positive results over time. The sensitivity and dilution management of strongly positive antibody levels were verified by titration of NIBSC 20/162.

Summary

By the correlation study TestLine confirmed high level of agreement between results obtained by Virus Neutralization Test and TL EIA COVID-19 RBD IgG kit. The level of neutralizing antibodies corresponds with level of anti-RBD IgG antibodies and therefore EIA kits developed by TestLine can be used to determine if concentration of antibodies has the protective effect. Moreover by harmonization with International Standards the kits can be used for monitoring the antibody concentration after vaccination by the time.

Ordering information:

| <u>Cat. No.</u> | <u>Product</u> | <u>No. of wells</u> |
|-----------------|---------------------------|---------------------|
| CoRA96 | EIA COVID-19 RBD IgA | 96 |
| CoRG96 | EIA COVID-19 RBD IgG | 96 |
| CoRM96 | EIA COVID-19 RBD IgM | 96 |
| SK-CoRA96 | SmartEIA COVID-19 RBD IgA | 96 |
| SK-CoRG96 | SmartEIA COVID-19 RBD IgG | 96 |
| SK-CoRM96 | SmartEIA COVID-19 RBD IgM | 96 |

SmartEIA kits are designed for automated processing using the Agility® analyser.

TestLine Clinical Diagnostics Ltd.

Krizikova 68, 612 00 Brno, Czech Republic

+420 549 121 203

sales@testlinecd.com

www.testlinecd.com