

Gemini

Compact, fully automated ELISA and MBA analyzer



strattec●●


 **TestLine®**

Full automation of ELISA method

- Fully open, 3 plate system automatic analyzer
- Highly flexible time manager with LIS connectivity
- Easy installation and maintenance
- High pipetting safety ensured by triple control
- Convenient programming of test protocols
- Easy operation
- Proven reliability
- Possibility of online support
- The system complies with the diagnostic guidelines for in vitro diagnostics


Incubation

- Two independent incubators (to 45 °C) with shaking (20 Hz)
- Three incubators at RT




Photometer

- Spectral photometer with up to 8 filters
- Reading of one or two wavelengths
- End point and kinetics
- Shaking



Washing

- 8-channel washing station (aspiration and dosing)
- Three wash solutions
- Multipoint washing, wetting, top and bottom washing, different aspiration speed




Simplicity and flexibility

- Flexible scheduling of two plates
- Graphical representation of scheduled operations
- Proven software reliability
- LIS connectivity
- Established software platform installed in more than 9,000 locations worldwide

Safety and processes control

- Pipetting pressure control
- Detection of clots and control of sufficient sample volume
- Colorimetric verification of dosing accuracy
- Identification of samples and reagents barcodes





Capacity and throughput

- Initial loading of up to 192 samples
- Continuous sample loading and retrieval
- Three positions for tips and two positions for pre-dilution strips



Pipetting

- Liquid pipettor for disposable tips (300 µl or 1100 µl)
- Flexible programming of single or multi-step dilution, mixing or sample storage



Technical specification

Basic

Dimensions (W x H x D)	125 x 90 x 75 cm
-------------------------------	------------------

Pipettor

Min. / max. volumes	10 µl – 300 µl (300 µl tip) / 301 µl – 1000 µl (1100 µl tip)
----------------------------	-----------------------------------------------------------------

Dosing accuracy (one-step)	< 3% CV (20 µl) / < 3% CV (100 µl)
-----------------------------------	------------------------------------

Dosing accuracy (multi-step)	< 10% CV (16 x 20 µl) / < 3% CV (8 x 100 µl)
-------------------------------------	-------------------------------------------------

Highlights

pressure monitoring, capacitive detection of liquid levels, tips detection, mixing, multiple dilution steps, archiving

ELISA

Capacity	up to 192 samples flexibility: ex. 144 samples + 8 reagents + 16 controls
-----------------	---------------------------------------------------------------------------------

Photometer

Photometric range	0 – 3.0 OD
--------------------------	------------

Spectral range	400 – 700 nm (up to 8 filters)
-----------------------	--------------------------------

Reading time	< 15 sec, < 30 sec (2 wavelengths)
---------------------	------------------------------------

Range	1% CV (1.0 OD)
--------------	----------------

Linearity	0 – 2.000 OD ±1%
------------------	------------------

Incubator

Temperature range	to 45 °C
--------------------------	----------

Temperature distribution	±1,5 °C (temperature monitoring)
---------------------------------	----------------------------------

Shaking	20 Hz
----------------	-------

Washer

Accuracy	10% CV (300 µl)
-----------------	-----------------

Residual volume	< 2.5 µl U-plates
------------------------	-------------------

Washing capacity	3 washing solutions
-------------------------	---------------------

Highlights

multipoint washing, soaking, top and bottom washing, different aspiration speed

All values are achieved under optimum conditions and may vary depending on environmental, instrument and processing conditions.

MANUFACTURER:

strattec ● ●

Strattec SE
Gewerbestr. 37
75217 Birkenfeld
Germany

DISTRIBUTOR:

BI TestLine®

TestLine Clinical Diagnostics Ltd.
Krizikova 68, 612 00 Brno, Czech Republic
www.testlinecd.com

Sales Department:
+420 549 121 203
sales@testlinecd.com

Orders:
+420 541 243 390
order@testlinecd.com