

Name:	Anchor® Tips, 300 uL
Purpose:	Vessels used to transport sample from the laboratory container to the cuvette with an automated pipetting system in the KleeYa Evaluation Unit single use
Ident No.:	Anchor® Tips, 300 uL – Ident No.:10036506
Packaging units:	Anchor® Tips, 300 uL, Set of 5760 pieces
Material:	PP medical grade
BASIC UDI-DI:	91201239109889J4

1 INTENDED USE

The intended use is as a vessel for the transfer of liquids.

The Anchor® Tips are intended exclusively for single use only and for application by qualified staff.

Anchor® Tips are designed to be used with an automated pipetting system in the KleeYa instrument in line with *in vitro* diagnostic analyses.

INFO

The Anchor Tips are validated for use in combination with the KleeYa instrument. In their function of KleeYa accessories, they are classified as Class A IVD medical devices in regards to EU IVD regulation 2017/746. Conformity with the relevant General Safety and Performance Requirements is demonstrated for the Anchor Tips in connection with the KleeYa instrument.

For more information about the Anchor® Tips, 300 uL please refer to the KleeYa User Manual.

⚠ WARNING



Handling

This instruction is provided for your safety and gives important instructions for the handling of the article.

- Read all instructions (including KleeYa User Manual)!
- Do not use damaged articles!
- Do not use any articles whose direct packaging is damaged!
- Do not use any articles whose expiration date has been exceeded!
- Never reuse single-use articles! Repeatedly use of single-use articles might cause cross-contamination. Furthermore, correct product performance is not guaranteed anymore.

2 SPECIFICATIONS

Nominal volume: 300 µL

Nominal length: 55 mm

3 USE


With two independent drawers for Anchor® Tips trays, a continuous refilling of tips and therefore an uninterrupted operation of the KleeYa System is possible.


Every drawer can be filled with one Anchor® Tip tray with 96 tips each. The KleeYa System's software monitors the loading of the trays.

The Anchor® Tips trays are inserted into the Tip Loading Bay of the KleeYa System. The Tip Loading Bay has the capacity to hold up to two tip trays.

For further instructions for use see KleeYa User Manual "3.3.3 Load Or Unload Anchor® Tips".

4 ENVIRONMENTAL CONDITIONS

	Temperature	Operating: 18°C to 37°C (64,4°F to 98,6°F)
		Storage: 2°C to 32°C (35,6°F to 89,6°F) for a maximum of 3 years
		Transport: -10°C to 60°C (14°F to 140°F) for a maximum of 3 days

	Humidity	Operating: 20 – 80 % non-condensing
		Storage: 20 - 80 % non-condensing for a maximum of 3 years
		Transport: 10 - 90 % non-condensing for a maximum of 3 days

NOTICE Protect from sunlight and UV rays!

5 DISPOSAL AND DECONTAMINATION

NOTICE

Disposal of non-contaminated parts

Material out of use (e.g. packaging material, units with damaged direct packaging) should be properly disposed according to local waste disposal regulations.

⚠ WARNING

Disposal of contaminated parts

Potential infectious material will cause severe environmental contamination.

Strictly follow the local and national provisions, legislation and laboratory regulations.



PACKAGE INSERT

ANCHOR® TIPS, 300 µL

6 SUPPORT

NOTICEE-mail: Support.KleeYa@diatron.com

7 MANUFACTURER



STRATEC Consumables GmbH

Sonystrasse 20

5081 Anif

Austria

Internet: www.strattec.com/solutions/consumables