

## Dynex Agility System



**Catalog number:** 67000

**Storage:** Bez teplotního režimu

**Producer:** Dynex Technologies, Inc.

### Instrument description:

- Fully automated processing eliminating all manual liquid-transfer steps.
- Up to 16 SmartKits® carriers stored on-board for simultaneous runs.
- Flexible throughput allows up to 12 plates on-board at once – the instrument utilizes three robotic arms to obtain maximum process efficiency.
- Up to 200 sample tubes on-board with additional continuous load capability.
- Comprehensive consumables and inventory monitoring system enabling continual reagents handling.
- Automated kits' barcodes scanning, automatic tube type recognition and automatic strip check.
- Comprehensive LIS integration.
- Intuitive, easy-to-use interface.
- All consumables for each assay are already prepackaged in the TestLine SmartKit® for the maximum user comfort and safety.