

AMP Rapid Test SARS-CoV-2 Ag (10 cassettes)



Catalog number: RT2951

Storage: RT

Producer: AMEDA Labordiagnostik GmbH

- **AMP Rapid Test SARS-CoV-2 Ag** is a rapid immuno-chromatographic test for qualitative detection of SARS-CoV-2 nucleocapsid protein antigen in human nasal or throat swab samples as an aid in quick and efficient diagnosis of Corona Virus Disease 2019 (COVID-19)
- Clinical studies evidenced high sensitivity and specificity of the test in comparison with RT-PCR. Potential cross-reactivity was tested with a variety of common pathogens and the concentrations at which no cross-reactivity was observed are provided in the IFU.
- Kit size - 10 cassettes.
Sterile swabs for collection of nasopharyngeal samples, as well as buffer and tubes for sample extraction are included in the kits.
- **Reading time:** 15 to max. 20 minutes
Detection limit: 1.15×10^2 TCID₅₀/mL

Performance: compared with RT-PCR

Test Sensitivity:	97,3 %	(95% CI: 90.0% - 99.8%)
Test Specificity:	100,00 %	(95% CI: 96.6% - 100%)
Overall agreement:	98,8 %	(95% CI: 91.8% - 99.9%)

- No special equipment required.

Please contact individual area sales managers for your country for further information.