

Antibody Index Software

Catalog number:

Storage: Bez teplotního režimu

Antibody index – software enables evaluation of antibody index (AI), it is ratio of specific antibodies in cerebrospinal fluid and in serum and state of hematoliquor barrier and total immunoglobulin concentration in cerebrospinal fluid and serum, using diagnostic products TestLine.

Software information:

- Effective and user friendly solution for reliable guidance through the evaluation.
- Calibration curve incorporaed into the software.
- The software enables saving and printing of two Results types summary or individual.
- All results are archived in the software database.
- Enables connectivity to laboratory information system.

Current version of the software: 2.9.

Exportováno 26. 10. 2025 Stránka 1 / 1