

## Antibody Index Software

**Catalog number:**

**Storage:** Bez teplotního režimu

Antibody index – software enables evaluation of antibody index (AI), it is ratio of specific antibodies in cerebrospinal fluid and in serum and state of hemato-liquor barrier and total immunoglobulin concentration in cerebrospinal fluid and serum, using diagnostic products TestLine.

### Software information:

- Effective and user friendly solution for reliable guidance through the evaluation.
- Calibration curve incorporated into the software.
- The software enables saving and printing of two Results types - summary or individual.
- All results are archived in the software database.
- Enables connectivity to laboratory information system.

**Current version of the software: 2.9.**