

## Antibody Index Software

**EAN Code:** 8595635304996

**Catalog number:** SwAI01

**Package size:** 1 pcs

**Storage:** Bez teplotního režimu

**Producer:** TestLine Clinical Diagnostics s.r.o.

Antibody Index Software enables evaluation of the antibody index (AI), i.e. the ratio of specific antibodies in cerebrospinal fluid (CSF) and serum in relation to the status of the blood–CSF barrier and the concentration of total immunoglobulins in CSF and serum, using TestLine diagnostic products.

### Program information:

- User-friendly and intuitively guides the user through the entire evaluation process.
- Calibration curve integrated into the software.
- The program allows saving and printing of two types of evaluations – summary and individual.
- All results are archived in the program database.