

## **CLIA Borrelia CSF IgM**



EAN Code: 8595635312144 Catalog number: CL-BCSFM50

Package size: 50 tests

Storage: 2-8 °C

**Producer:** TestLine Clinical Diagnostics s.r.o.

The CLIA imunoassay is based on the principle of magnetic particle separation, and the detection system utilizes flash chemiluminescence with the help of acridinium ester, and the intensity of the light flash is measured in relative light units (RLU).

## General information about the kit:

- Each individual test contains all the necessary reagents (ready-to-use) for conducting the examination.
- Calibration is included in the package.
- Results are evaluated using U/ml.

For more information, please refer to **BiovendorGroup CLIA** webpage.

Exportováno 26. 10. 2025 Stránka 1 / 1