

CLIA CRP wide-range

EAN Code: 8595635311413
Catalog number: CL-CRP100
Package size: 100 tests
Storage: 2-8 °C
Producer: TestLine Clinical Diagnostics s.r.o.



The CLIA immunoassay is based on the principle of magnetic particle separation, and the detection system utilizes flash chemiluminescence with the help of acridinium ester, and the intensity of the light flash is measured in relative light units (RLU).

General information about the kit:

- Each individual test contains all the necessary reagents (ready-to-use) for conducting the examination.
- Calibration is included in the package.
- Results are evaluated using mg/l.

For more information, please refer to [BiovendorGroup CLIA webpage](#).