

## **CLIA Tetanus Toxoid IgG**



EAN Code: 8595635311925 Catalog number: CL-TeTG100

Package size: 100 Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.

The CLIA imunoassay is based on the principle of magnetic particle separation, and the detection system utilizes flash chemiluminescence with the help of acridinium ester, and the intensity of the light flash is measured in relative light units (RLU).

## General information about the kit:

- Each individual test contains all the necessary reagents (ready-to-use) for conducting the examination.
- Calibration is included in the package.
- Results are evaluated using IU/ml.

For more information, please refer to BiovendorGroup CLIA webpage.

Exportováno 26. 10. 2025 Stránka 1 / 1