

Control set CLIA COVID-19 NP IgG

EAN Code: 8595635310782

Catalog number: CL-CoNGCON

Package size: 2x20 tests

Storage: 2-8 °C

Producer: TestLine Clinical Diagnostics s.r.o.



- The kit is intended for use together with CLIA kits on the KleeYa® analyzer.
- The kit includes positive and negative control sera, which are used to verify the functionality of CLIA tests and the accuracy of obtained results.

For more information, please refer to [BiovendorGroup CLIA webpage](#).