


## RoboBlot



**Catalog number:** RBT-001  
**Storage:** Bez teplotního režimu  
**Producer:** Bee Robotics Ltd.

- Fully automated performance of the assays (sample pipetting from primary tube, reagent dosing, incubation, strip washing and result evaluation directly in incubation tray)
- Full compatibility of RoboBlot, TestLine Immunoblot Software and TestLine kits
- Control of the entire process with Immunoblot Software and analyser's touch screen
- Connectivity to LIS
- Samples' barcode reader
- Integrated camera for strips scanning
- Minimal maintenance operations and hands on time necessary
- Capacity – 50 samples
- Flexibility – 5 different methods at the same time