TMB-Complete 2 (1 l)



EAN Code: 8595635303630
Catalog number: TMB0L1
Package size: 1000 ml
Storage: 2-8 °C
Producer: TestLine Clinical Diagnostics s.r.o.

Description:

- Solution contains TMB (3.3', 5.5'-TETRAMETHYLBENZIDINE), hydrogen peroxide (H_2O_2) and stabilizers.
- Solution does not contain any toxic compounds.
- The guarantied absorbance for 100 µl blank is less than 0.100 at 450 nm.
- Shelf life: 30 months.